

PRODUCT SPECIFICATION

Name β 2-Microglobulin Monoclonal Antibody (β 2-MG mAb)
Code HM059A

RECEIVER INFORMATION

Storage Temperature 2-8 °C (up to 2 Years), below -20 °C (up to 3 Years)
Storage Notes Avoid repeat freeze thaw cycles

PRODUCT INFORMATION

Description Monoclonal antibody directed to β 2-MG

Specificity Specific for β 2-MG

Information

Clone	Host /Source	Isotype	Purification
1059E	Mouse	IgG1	Protein G
1059F	Mouse	IgG1	Protein G
1059G	Mouse	IgG1	Protein G
1059H	Mouse	IgG1	Protein G

Application

Pair	Capture	Detection	Recommended use
Pair 1	1059E	1059F	CLIA/POCT
Pair 2	1059G	1059H	CLIA/POCT

Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Presentation Single homogenous batch, 20mM PBS, pH7.2, 0.05% Proclin300

HEALTH AND SAFETY

Warning For in-vitro use and research or further processing only.
This product contains Proclin300, which has been classified as Xi (Irritant), in European Directive 1999/45/EC. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water.