

## PRODUCT SPECIFICATION

**Name** 2019-nCoV N Protein Monoclonal Antibody  
(2019-nCoV NP mAb)

**Code** HM190BR

### RECEIVER INFORMATION

---

**Storage Temperature** 2-8 °C (up to 2 Years), below -20 °C (up to 3 Years), -80 °C (up to 5 Years)

**Storage Notes** Avoid repeat freeze thaw cycles

### PRODUCT INFORMATION

---

**Description** Monoclonal antibody directed to 2019-nCoV NP

**Specificity** Specific for 2019-nCoV NP

**Information**

Clone	Host /Source	Isotype	Purification
1190D	Rabbit	IgG	Protein G

**Application**

Pair	Capture	Detection	Recommended use
Pair 1	1190B	1190D	CLIA/POCT
Pair 2	1190D	1190B	CLIA/POCT

Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

**Presentation** Single homogenous batch, 20mM PBS, pH7.2, 0.05% Proclin300

### HEALTH AND SAFETY

---

**Warning**

For in-vitro use and research or further processing only.  
This product contains Proclin300, which has been classified as Xi (Irritant), in European Directive 1999/45/EC. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water.