

## PRODUCT SPECIFICATION

**Name** Adiponectin Monoclonal Antibody (ADPN mAb)  
**Code** HM073C

### RECEIVER INFORMATION

---

**Storage Temperature** 2-8 °C (up to 1 Month), below -20 °C (up to 18 Months)  
**Storage Notes** Avoid repeat freeze thaw cycles

### PRODUCT INFORMATION

---

**Description** Monoclonal antibody directed to ADPN

**Specificity** Specific for ADPN

**Information**

Clone	Host /Source	Isotype	Purification
1073D	Mouse	IgG	Affinity Chromatography

**Application**

Pair	Capture	Detection	Recommended use
Pair 1		1073D	Biochemical

Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

**Presentation** Single homogenous batch, 60nM Tris, pH8.1

### HEALTH AND SAFETY

---

**Warning** For in-vitro use and research or further processing only.