

PRODUCT SPECIFICATION

Name Carbohydrate Antigen 15-3 Monoclonal Antibody (CA15-3 mAb)

Code HM105B

RECEIVER INFORMATION

Storage Temperature 2-8 °C (up to 2 Years), below -20 °C (up to 3 Years), -80 °C (up to 5 Years)

Storage Notes Avoid repeat freeze thaw cycles

PRODUCT INFORMATION

Description Monoclonal antibody directed to CA15-3

Specificity Specific for CA15-3

Information

Clone	Host /Source	Isotype	Purification
1105D	HEK293	IgG1	Protein G
1105E	Mouse	IgG1	Protein G

Application

Pair	Capture	Detection	Recommended use
Pair 1	1105D	1105E	CLIA

Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Presentation Single homogenous batch, 20mM PBS, pH7.2, 0.05% Proclin300

HEALTH AND SAFETY

Warning

For in-vitro use and research or further processing only.
This product contains Proclin300, which has been classified as Xi (Irritant), in European Directive 1999/45/EC. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water.