

PRODUCT SPECIFICATION

Name Carbohydrate Antigen 72-4 Monoclonal Antibody (CA72-4 mAb)

Code HM107B

RECEIVER INFORMATION

Storage Temperature 2-8 °C (up to 2 Years), below -20 °C (up to 3 Years), -80 °C (up to 5 Years)

Storage Notes Avoid repeat freeze thaw cycles

PRODUCT INFORMATION

Description Monoclonal antibody directed to CA72-4

Specificity Specific for CA72-4

Information

Clone	Host /Source	Isotype	Purification
1107D	Mouse	IgG1	Protein G
1107E	Mouse	IgG1	Protein G

Application

Pair	Capture	Detection	Recommended use
Pair 1	1107D	1107E	CLIA

Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Presentation Single homogenous batch, 20mM PBS, pH7.2, 0.05% Proclin300

HEALTH AND SAFETY

Warning

For in-vitro use and research or further processing only.
This product contains Proclin300, which has been classified as Xi (Irritant), in European Directive 1999/45/EC. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water.