

## PRODUCT SPECIFICATION

<b>Name</b>	<b>Collagen Type IV Protein Monoclonal Antibody (CIV mAb)</b>
<b>Code</b>	<b>HM026B</b>

---

### RECEIVER INFORMATION

<b>Storage Temperature</b>	2-8 °C (up to 1 Month), below -20 °C (up to 5 Years)
<b>Storage Notes</b>	Avoid repeat freeze thaw cycles

---

### PRODUCT INFORMATION

<b>Description</b>	Monoclonal antibody directed to CIV
<b>Specificity</b>	Specific for CIV
<b>Information</b>	

Clone	Host /Source	Isotype	Purification
1026A	Mouse	IgG1	Protein A
1026B	Mouse	IgG1	Protein A

**Application**

Pair	Capture	Detection	Recommended use
Pair 1	1026A	1026B	CLIA/Biochemical

Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

<b>Presentation</b>	Single homogenous batch, 10mM PBS, pH7.4
---------------------	--

---

### HEALTH AND SAFETY

<b>Warning</b>	For in-vitro use and research or further processing only.
----------------	---