

PRODUCT SPECIFICATION

Name Cytokeratin 19 Fragment 21-1 Monoclonal Antibody (CYFRA21-1 mAb)

Code HM102A

RECEIVER INFORMATION

Storage Temperature 2-8 °C (up to 3 Years)

Storage Notes Avoid freezing

PRODUCT INFORMATION

Description Monoclonal antibody directed to CYFRA21-1

Specificity Specific for CYFRA21-1

Information

Clone	Host /Source	Isotype	Purification
1102A	Murine	IgG1	Ion exchange
1102B	Murine	IgG1	Ion exchange

Application

Pair	Capture	Detection	Recommended use
Pair 1	1102A	1102B	CLIA

Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Presentation Single homogenous batch, PBS, 0.05% Sodium Azide

HEALTH AND SAFETY

Warning

For in-vitro use and research or further processing only.
This product contains sodium azide, which has been classified as Xn (Harmful), in European Directive 67/548/EEC in the concentration range of 0.1-1.0%. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent azide build-up in drains.