

## PRODUCT SPECIFICATION

**Name** Ferritin Monoclonal Antibody (Fer mAb)  
**Code** HM092B

### RECEIVER INFORMATION

---

**Storage Temperature** 2-8 °C (up to 2 Years), below -20 °C (up to 3 Years)  
**Storage Notes** Avoid repeat freeze thaw cycles

### PRODUCT INFORMATION

---

**Description** Monoclonal antibody directed to Fer  
**Specificity** Specific for Fer  
**Information**

Clone	Host /Source	Isotype	Purification
1092A	Mouse	IgG2b	Protein G
1092B	Mouse	IgG2b	Protein G
1092C	Mouse	IgG2a	Protein G
1092E	Mouse	IgG2b	Protein G
1092F	Mouse	IgG2a	Protein G

**Application**

Pair	Capture	Detection	Recommended use
Pair 1	1092B	1092C	POCT
Pair 2		1092A	Biochemical
Pair 3	1092E	1092F	CLIA

Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

**Presentation** Single homogenous batch, 20mM PBS, pH7.2, 0.05% Proclin300

### HEALTH AND SAFETY

---

**Warning** For in-vitro use and research or further processing only.  
This product contains Proclin300, which has been classified as Xi (Irritant), in European Directive 1999/45/EC. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water.