

PRODUCT SPECIFICATION

Name Influenza A Nucleoprotein Monoclonal Antibody (FluA NP mAb)

Code HM185B

RECEIVER INFORMATION

Storage Temperature 2-8 °C (up to 3 Months), below -20 °C (up to 1 Year), -70 °C (up to 3 Years)

Storage Notes Avoid repeat freeze thaw cycles

PRODUCT INFORMATION

Description Monoclonal antibody directed to FluA NP

Specificity Specific for FluA NP

Information

Clone	Host /Source	Isotype	Purification
1185C	Mouse Ascites	IgG1	Ammonium sulfate precipitation
1185D	Mouse Ascites	IgG2a	Ammonium sulfate precipitation

Application

Pair	Capture	Detection	Recommended use
Pair 1	1185C	1185D	POCT

Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Presentation Single homogenous batch, PBS, pH7.4, 0.1% Sodium Azide

HEALTH AND SAFETY

Warning For in-vitro use and research or further processing only. This product contains sodium azide, which has been classified as Xn (Harmful), in European Directive 67/548/EEC in the concentration range of 0.1-1.0%. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent azide build-up in drains.