

## PRODUCT SPECIFICATION

**Name** Influenza B Nucleoprotein Monoclonal Antibody (FluB NP mAb)

**Code** HM193B

### RECEIVER INFORMATION

---

**Storage Temperature** 2-8 °C (up to 3 Months), below -20 °C (up to 1 Year), -70 °C (up to 3 Years)

**Storage Notes** Avoid repeat freeze thaw cycles

### PRODUCT INFORMATION

---

**Description** Monoclonal antibody directed to FluB NP

**Specificity** Specific for FluB NP

**Information**

Clone	Host /Source	Isotype	Purification
1193C	Mouse Ascites	IgM	Ammonium sulfate precipitation
1193D	Mouse Ascites	IgG2b	Ammonium sulfate precipitation

**Application**

Pair	Capture	Detection	Recommended use
Pair 1	1193C	1193D	POCT

Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

**Presentation** Single homogenous batch, PBS, pH7.4, 0.1% Sodium Azide

### HEALTH AND SAFETY

---

**Warning** For in-vitro use and research or further processing only. This product contains sodium azide, which has been classified as Xn ( Harmful ), in European Directive 67/548/EEC in the concentration range of 0.1-1.0%. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent azide build-up in drains.