

PRODUCT SPECIFICATION

Name Hemoglobin A1c Monoclonal Antibody (HbA1c mAb)
Code HM004B

RECEIVER INFORMATION

Storage Temperature 2-8 °C (up to 2 Weeks), below -20 °C (up to 2 Years)
Storage Notes Avoid repeat freeze thaw cycles

PRODUCT INFORMATION

Description Monoclonal antibody directed to HbA1c

Specificity Specific for HbA1c

Information

Clone	Host /Source	Isotype	Purification
1004B	Mouse	IgG1	Protein G

Application

Pair	Capture	Detection	Recommended use
Pair 1	1004B	1016B	POCT

Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Presentation Single homogenous batch, 10mM PBS, pH7.4, 0.1% Sodium Azide

HEALTH AND SAFETY

Warning For in-vitro use and research or further processing only.
This product contains sodium azide, which has been classified as Xn (Harmful), in European Directive 67/548/EEC in the concentration range of 0.1-1.0%. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent azide build-up in drains.