

## PRODUCT SPECIFICATION

**Name** Hemoglobin A1c Monoclonal Antibody (HbA1c mAb)  
**Code** HM004E

### RECEIVER INFORMATION

---

**Storage Temperature** 2-8 °C (up to 2 Years), below -20 °C (up to 3 Years), -80 °C (up to 5 Years)  
**Storage Notes** Avoid repeat freeze thaw cycles

### PRODUCT INFORMATION

---

**Description** Monoclonal antibody directed to HbA1c  
**Specificity** Specific for HbA1c. No cross-reactivity with human Hb  
**Information**

Clone	Host /Source	Isotype	Purification
1004F	Mouse	IgG1	Protein G
1004G	Mouse	IgG2a	Protein G

**Application**

Pair	Capture	Detection	Recommended use
Pair 1		1004F	Biochemical
Pair 2		1004G	Biochemical

Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

**Presentation** Single homogenous batch, 20mM PBS, pH7.2, 0.05% Proclin300

### HEALTH AND SAFETY

---

**Warning** For in-vitro use and research or further processing only.  
This product contains Proclin300, which has been classified as Xi (Irritant), in European Directive 1999/45/EC. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water.