

PRODUCT SPECIFICATION

Name Interferon-gamma Monoclonal Antibody (IFN-gamma mAb)

Code HM211A

RECEIVER INFORMATION

Storage Temperature 2-8 °C (up to 3 Months), below -20 °C (up to 1 Year), -70 °C (up to 3 Years)

Storage Notes Avoid repeat freeze thaw cycles

PRODUCT INFORMATION

Description Monoclonal antibody directed to IFN-gamma

Specificity Specific for IFN-gamma

Information

Clone	Host /Source	Isotype	Purification
1211A	Mouse Ascites	IgG1	Ammonium sulfate precipitation
1211B	Mouse Ascites	IgG1	Ammonium sulfate precipitation

Application

Pair	Capture	Detection	Recommended use
Pair 1	1211A	1211B	ELISA
Pair 2	1211B	1211A	ELISA

Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Presentation Single homogenous batch, Saline, 0.1% Sodium Azide

HEALTH AND SAFETY

Warning

For in-vitro use and research or further processing only. This product contains sodium azide, which has been classified as Xn (Harmful), in European Directive 67/548/EEC in the concentration range of 0.1-1.0%. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent azide build-up in drains.