

PRODUCT SPECIFICATION

Name Insulin Like Growth Factor Binding Protein-1
Monoclonal Antibody (IGFBP-1 mAb)

Code HM041A

RECEIVER INFORMATION

Storage Temperature 2-8 ° C (up to 2 Years), below -20 ° C (up to 3 Years)

Storage Notes Avoid repeat freeze thaw cycles

PRODUCT INFORMATION

Description Monoclonal antibody directed to IGFBP-1

Specificity Specific for IGFBP-1

Information

Clone	Host /Source	Isotype	Purification
1041A	Mouse	IgG1	Protein G
1041B	Mouse	IgG2b	Protein G
1041C	Mouse	IgG2b	Protein G

Application

Pair	Capture	Detection	Recommended use
Pair 1	1041A	1041B	POCT
Pair 2	1041A	1041C	POCT

Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Presentation Single homogenous batch, 20mM PBS, pH7.2, 0.05% Proclin300

HEALTH AND SAFETY

Warning For in-vitro use and research or further processing only.
This product contains Proclin300, which has been classified as Xi (Irritant), in European Directive 1999/45/EC. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water.