

## PRODUCT SPECIFICATION

**Name** Myoglobin Monoclonal Antibody (Mb mAb)  
**Code** HM035A

### RECEIVER INFORMATION

---

**Storage Temperature** 2-8 °C (up to 2 Years), below -20 °C (up to 3 Years), -80 °C (up to 5 Years)  
**Storage Notes** Avoid repeat freeze thaw cycles

### PRODUCT INFORMATION

---

**Description** Monoclonal antibody directed to Mb  
**Specificity** Specific for Mb  
**Information**

Clone	Host /Source	Isotype	Purification
1035C	Mouse	IgG2b	Protein G
1035D	Mouse	IgG1	Protein G
1035E	Mouse	IgG1	Protein G
1035F	Mouse	IgG1	Protein G

**Application**

Pair	Capture	Detection	Recommended use
Pair 1	1035C	1035D	POCT
Pair 2	1035E	1035F	CLIA

Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

**Presentation** Single homogenous batch, 20mM PBS, pH7.2, 0.05% Proclin300

### HEALTH AND SAFETY

---

**Warning** For in-vitro use and research or further processing only.  
This product contains Proclin300, which has been classified as Xi (Irritant), in European Directive 1999/45/EC. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water.