

PRODUCT SPECIFICATION

Name Pepsinogen I Monoclonal Antibody (PGI mAb)
Code HM023A

RECEIVER INFORMATION

Storage Temperature 2-8 °C (up to 2 Years), below -20 °C (up to 3 Years)
Storage Notes Avoid repeat freeze thaw cycles

PRODUCT INFORMATION

Description Monoclonal antibody directed to PGI
Specificity Specific for PGI
Information

Clone	Host /Source	Isotype	Purification
1023A	Mouse	IgG1	Protein G
1023B	Mouse	IgG1	Protein G
1023E	Mouse	IgG1	Protein G
1023F	Mouse	IgG1	Protein G
1023G	Mouse	IgG2b	Protein G

Application

Pair	Capture	Detection	Recommended use
Pair 1	1023A	1023B	Biochemical
Pair 2	1023E	1023F	POCT
Pair 3	1023F	1023G	CLIA

Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Presentation Single homogenous batch, 20mM PBS, pH7.2, 0.05% Proclin 300

HEALTH AND SAFETY

Warning For in-vitro use and research or further processing only.
This product contains Proclin300, which has been classified as Xi (Irritant), in European Directive 1999/45/EC. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water.