

PRODUCT SPECIFICATION

Name Pepsinogen II Monoclonal Antibody (PGII mAb)
Code HM024A

RECEIVER INFORMATION

Storage Temperature 2-8 °C (up to 2 Years), below -20 °C (up to 3 Years), -80 °C (up to 5 Years)
Storage Notes Avoid repeat freeze thaw cycles

PRODUCT INFORMATION

Description Monoclonal antibody directed to PGII

Specificity Specific for PGII

Information

Clone	Host /Source	Isotype	Purification
1024A	Mouse	IgG1	Protein G
1024B	Mouse	IgG1	Protein G
1024E	Mouse	IgG1	Protein G
1024F	Mouse	IgG1	Protein G
1024G	Mouse	IgG2b	Protein G

Application

Pair	Capture	Detection	Recommended use
Pair 1	1024A	1024B	Biochemical
Pair 2	1024E	1024F	POCT
Pair 3	1024A	1024G	Biochemical/CLIA

Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Presentation Single homogenous batch, 20mM PBS, pH7.2, 0.05% Proclin 300

HEALTH AND SAFETY

Warning

For in-vitro use and research or further processing only.
This product contains Proclin300, which has been classified as Xi (Irritant), in European Directive 1999/45/EC. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water