

## PRODUCT SPECIFICATION

**Name** Serum Amyloid A Protein Monoclonal Antibody (SAA mAb)

**Code** HM021I

### RECEIVER INFORMATION

---

**Storage Temperature** 2-8 °C (up to 2 Years), below -20 °C (up to 3 Years)

**Storage Notes** Avoid repeat freeze thaw cycles

### PRODUCT INFORMATION

---

**Description** Monoclonal antibody directed to SAA

**Specificity** Specific for SAA

**Information**

Clone	Host /Source	Isotype	Purification
1021D	Mouse	IgG2a	Protein G
1021E	Mouse	IgG1	Protein G
1021F	Mouse	IgG1	Protein G

**Application**

Pair	Capture	Detection	Recommended use
Pair 1	1021E	1021F	Biochemical
Pair 2	1021D	1021E	POCT

Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

**Presentation** Single homogenous batch, 20mM PBS, pH7.2, 0.05% Proclin300

### HEALTH AND SAFETY

---

**Warning**

For in-vitro use and research or further processing only.  
This product contains Proclin300, which has been classified as Xi (Irritant), in European Directive 1999/45/EC. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water.