

## PRODUCT SPECIFICATION

**Name** Triiodothyronin Monoclonal Antibody (T3 mAb)  
**Code** HM150A

### RECEIVER INFORMATION

---

**Storage Temperature** 2-8 °C (up to 2 Years), below -20 °C (up to 3 Years)  
**Storage Notes** Avoid repeat freeze thaw cycles

### PRODUCT INFORMATION

---

**Description** Monoclonal antibody directed to T3

**Specificity** Specific for T3

**Information**

| Clone | Host /Source | Isotype | Purification |
|-------|--------------|---------|--------------|
| 1150A | Goat         | IgG1    | Protein G    |

**Application**

| Pair   | Capture | Detection | Recommended use |
|--------|---------|-----------|-----------------|
| Pair 1 |         | 1150A     | CLIA            |

Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

**Presentation** Single homogenous batch, 20mM PBS, pH7.2, 0.05% Proclin300

### HEALTH AND SAFETY

---

**Warning** For in-vitro use and research or further processing only.  
This product contains Proclin300, which has been classified as Xi (Irritant), in European Directive 1999/45/EC. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water.