

PRODUCT SPECIFICATION

Name Thyroid Peroxidase Monoclonal Antibody (TPO mAb)
Code HM187A

RECEIVER INFORMATION

Storage Temperature 2-8 °C (up to 2 Years), below -20 °C (up to 3 Years)
Storage Notes Avoid repeat freeze thaw cycles

PRODUCT INFORMATION

Description Monoclonal antibody directed to TPO

Specificity Specific for TPO

Information

Clone	Host /Source	Isotype	Purification
1187A	Mouse	IgG1	Protein G
1187B	Mouse	IgG2a	Protein G
1187C	Mouse	IgG1	Protein G
1187D	Mouse	IgG2a	Protein G

Application

Pair	Capture	Detection	Recommended use
Pair 1	1187A	1187B	ELISA
Pair 1	1187C	1187D	ELISA

Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Presentation Single homogenous batch, 20mM PBS, pH7.2, 0.05% Proclin300

HEALTH AND SAFETY

Warning For in-vitro use and research or further processing only.
This product contains Proclin300, which has been classified as Xi (Irritant), in European Directive 1999/45/EC. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water.