

PRODUCT SPECIFICATION

Name Free Prostate Special Antigen Monoclonal Antibody (free-PSA mAb)

Code HM168A

RECEIVER INFORMATION

Storage Temperature 2-8 °C (up to 3 Months), below -20 °C (up to 1 Year), -70 °C (up to 3 Years)

Storage Notes Avoid repeat freeze thaw cycles

PRODUCT INFORMATION

Description Monoclonal antibody directed to free-PSA

Specificity Specific for free-PSA

Information

Clone	Host /Source	Isotype	Purification
1168A	Mouse Ascites	IgG1	Protein A
1168B	Mouse Ascites	IgG1	Protein A

Application

Pair	Capture	Detection	Recommended use
Pair 1	1168A	1168B	CLIA/POCT

Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Presentation Single homogenous batch, PBS, pH7.2-7.4, 0.09% Sodium Azide

HEALTH AND SAFETY

Warning

For in-vitro use and research or further processing only.
This product contains sodium azide, which has been classified as Xn (Harmful), in European Directive 67/548/EEC in the concentration range of 0.1-1.0%. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent azide build-up in drains.