

## PRODUCT SPECIFICATION

**Name** Heart Type Fatty Acid Binding Protein Monoclonal Antibody (hFABP mAb)

**Code** HM040B

### RECEIVER INFORMATION

---

**Storage Temperature** 2-8 °C (up to 2 Years), below -20 °C (up to 3 Years)

**Storage Notes** Avoid repeat freeze thaw cycles

### PRODUCT INFORMATION

---

**Description** Monoclonal antibody directed to hFABP

**Specificity** Specific for hFABP

**Information**

Clone	Host /Source	Isotype	Purification
1040C	Mouse	IgG1	Protein G
1040D	Mouse	IgG1	Protein G
1040E	Mouse	IgG1	Protein G
1040F	Mouse	IgG2b	Protein G

**Application**

Pair	Capture	Detection	Recommended use
Pair 1	1040C	1040D	Biochemical
Pair 2	1040E	1040F	POCT

Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

**Presentation** Single homogenous batch, 20mM PBS, pH7.2, 0.05% Proclin300

### HEALTH AND SAFETY

---

**Warning** For in-vitro use and research or further processing only.  
This product contains Proclin300, which has been classified as Xi (Irritant), in European Directive 1999/45/EC. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water.