

## PRODUCT SPECIFICATION

<b>Name</b>	<b>Apolipoprotein A1 Antigen (ApoA1 Ag)</b>
<b>Code</b>	<b>HA075AN</b>

---

### RECEIVER INFORMATION

<b>Storage Temperature</b>	2-8 °C (up to 6 Months), below -20 °C (up to 2 Years)
<b>Storage Notes</b>	Avoid repeat freeze thaw cycles

---

### PRODUCT INFORMATION

<b>Description</b>	ApoA1 Ag
<b>Specificity</b>	Specific for ApoA1 antibody
<b>Application</b>	Suitable for use in Immunogen/Calibrators/Control. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
<b>Host /Source</b>	Native
<b>Presentation</b>	Single homogenous batch, 20mM PBS, pH7.2, 0.05% Proclin300

---

### HEALTH AND SAFETY

<b>Warning</b>	For in-vitro use and research or further processing only. This product contains Proclin300, which has been classified as Xi (Irritant), in European Directive 1999/45/EC. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water.
----------------	---