

PRODUCT SPECIFICATION

Name	Apolipoprotein B Antigen (ApoB Ag)
Code	HA076AN

RECEIVER INFORMATION

Storage Temperature	2-8 °C (up to 1 Year)
Storage Notes	Avoid freezeing

PRODUCT INFORMATION

Description	ApoB Ag
Specificity	Specific for ApoB antibody
Application	Suitable for use in Immunogen/Calibrators/Control. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Host /Source	Native
Presentation	Single homogenous batch, 20mM PBS, pH7.2, 0.05% Proclin300

HEALTH AND SAFETY

Warning	For in-vitro use and research or further processing only. This product contains Proclin300, which has been classified as Xi (Irritant), in European Directive 1999/45/EC. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water.
----------------	---