

PRODUCT SPECIFICATION

Name	Complement C3 Antigen (C3 Ag)
Code	HA087AN

RECEIVER INFORMATION

Storage Temperature	2-8 °C (up to 1 Year), below -20 °C (up to 2 Years)
Storage Notes	Avoid repeat freeze thaw cycles

PRODUCT INFORMATION

Description	C3 Ag
Specificity	Specific for C3 antibody
Application	Suitable for use in Immunogen/Calibrators/Control. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Host /Source	Native
Presentation	Single homogenous batch, 20mM PB, 1M NaCl, pH7.4, 0.05% Proclin300

HEALTH AND SAFETY

Warning	For in-vitro use and research or further processing only. This product contains Proclin300, which has been classified as Xi (Irritant), in European Directive 1999/45/EC. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water.
----------------	---