

PRODUCT SPECIFICATION

Name	Carbohydrate Antigen 72-4 Antigen (CA72-4 Ag)
Code	HA107AN

RECEIVER INFORMATION

Storage Temperature	Below -20 °C (up to 3 Years)
Storage Notes	Avoid repeat freeze thaw cycles

PRODUCT INFORMATION

Description	CA72-4 Ag
Specificity	Specific for CA72-4 antibody
Application	Suitable for use in Immunogen/Calibrators. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Host /Source	Native
Presentation	Single homogenous batch, 100mM PBS, pH7.2-7.4, 3% Sucrose, 0.1% Sodium Azide

HEALTH AND SAFETY

Warning	For in-vitro use and research or further processing only. This product contains sodium azide, which has been classified as Xn (Harmful), in European Directive 67/548/EEC in the concentration range of 0.1-1.0%. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent azide build-up in drains.
----------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------