

PRODUCT SPECIFICATION

| | |
|-------------|--|
| Name | Creatine Kinase Isoenzyme MB Antigen (CK-MB Ag) |
| Code | HA045BN |

RECEIVER INFORMATION

| | |
|----------------------------|--|
| Storage Temperature | 2-8 °C (up to 2 Weeks), below -20 °C (up to 2 Years) |
| Storage Notes | Avoid repeat freeze thaw cycles |

PRODUCT INFORMATION

| | |
|---------------------|---|
| Description | CK-MB Ag |
| Specificity | Specific for CK-MB antibody |
| Application | Suitable for use in Calibrators/Control. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded. |
| Host /Source | Native |
| Presentation | Single homogenous batch, Tris-Hcl, pH7.2 |

HEALTH AND SAFETY

| | |
|----------------|---|
| Warning | For in-vitro use and research or further processing only. This product contains Proclin300, which has been classified as Xi (Irritant), in European Directive 1999/45/EC. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water. |
|----------------|---|