

## PRODUCT SPECIFICATION

**Name** Cytokeratin 19 Fragment 21-1 Antigen  
(CYFRA21-1 Ag)

**Code** HA102AN

---

### RECEIVER INFORMATION

**Storage Temperature** Below -20 °C (up to 3 Years)

**Storage Notes** Avoid repeat freeze thaw cycles

---

### PRODUCT INFORMATION

**Description** CYFRA21-1 Ag

**Specificity** Specific for CYFRA21-1 antibody

**Application** Suitable for use in Immunogen.  
Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

**Host /Source** Native

**Presentation** Single homogenous batch, 100mM PBS, 0.05% Sodium Azide

---

### HEALTH AND SAFETY

**Warning** For in-vitro use and research or further processing only.  
This product contains sodium azide, which has been classified as Xn ( Harmful ), in European Directive 67/548/EEC in the concentration range of 0.1-1.0%. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent azide build-up in drains.