

PRODUCT SPECIFICATION

Name	Immunoglobulin G Antigen (IgG Ag)
Code	HA085DN

RECEIVER INFORMATION

Storage Temperature	2-8 °C (up to 1 Year), below -20 °C (up to 2 Years)
Storage Notes	Avoid repeat freeze thaw cycles

PRODUCT INFORMATION

Description	IgG Ag
Specificity	Specific for IgG antibody
Application	Suitable for use in Calibrators/Control. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Host /Source	Native
Presentation	Single homogenous batch, PBS, pH7.4, 0.05% Proclin 300

HEALTH AND SAFETY

Warning	For in-vitro use and research or further processing only. This product contains Proclin300, which has been classified as Xi (Irritant), in European Directive 1999/45/EC. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water.
----------------	---