

PRODUCT SPECIFICATION

Name Myeloperoxidase Antigen (MPO Ag)
Code HA039BN

RECEIVER INFORMATION

Storage Temperature Below -20 ° C (up to 3 Years)
Storage Notes Avoid repeat freeze thaw cycles

PRODUCT INFORMATION

Description MPO Ag
Specificity Specific for MPO antibody
Application Suitable for use in Immunogen/Calibrators
Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Host /Source Native
Presentation Single homogenous batch, 20mM Tris, 100mM NaCl, pH7.5

HEALTH AND SAFETY

Warning For in-vitro use and research or further processing only.