

PRODUCT SPECIFICATION

Name	Pepsinogen I Antigen (PGI Ag)
Code	HA023CN

RECEIVER INFORMATION

Storage Temperature	Below -20 °C (up to 3 Years)
Storage Notes	Avoid repeat freeze thaw cycles

PRODUCT INFORMATION

Description	PGI Ag
Specificity	Specific for PGI antibody
Application	Suitable for use in Calibrators/Control. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Host /Source	Human gastric mucosa
Presentation	Single homogenous batch, Tris-HCl, pH8.5

HEALTH AND SAFETY

Warning	For in-vitro use and research or further processing only.
----------------	---