

PRODUCT SPECIFICATION

Name	Rheumatoid Factor Antigen (RF Ag)
Code	HA129AN

RECEIVER INFORMATION

Storage Temperature	2-8 °C (up to 1 Week), below -20 °C (up to 1 Year)
Storage Notes	Avoid repeat freeze thaw cycles

PRODUCT INFORMATION

Description	RF Ag
Specificity	Specific for RF antibody
Application	Suitable for use in POCT/LETIA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Host /Source	Native
Presentation	Single homogenous batch, 20mM PBS, pH7.4

HEALTH AND SAFETY

Warning	For in-vitro use and research or further processing only.
----------------	---