

PRODUCT SPECIFICATION

Name	Cardiac Troponin I Antigen (cTnI Ag)
Code	HA033AN

RECEIVER INFORMATION

Storage Temperature	2-8 °C (up to 2 Weeks), below -20 °C (up to 3 Years)
Storage Notes	Avoid repeat freeze thaw cycles

PRODUCT INFORMATION

Description	cTnI Ag
Specificity	Specific for cTnI antibody
Application	Suitable for use in Calibrators/Control/CLIA/LETIA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Host /Source	Native
Presentation	Single homogenous batch, Tris, 6M Urea

HEALTH AND SAFETY

Warning	For in-vitro use and research or further processing only.
----------------	---